



Nalgene® square PETG media bottles

Product Images

Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 250 mL Capacity (anglo): 8 oz. Diameter (anglo): 2.4 in. Inner diameter (Anglo-American) neck: 1.10 in. For use with (application): Sampling, Preparation, Containment, Prolonged Storage Height (metric) with closure: 144 mm Number per carton: 48 Number per package: 12 Temperature range (metric): gt;40deg;C Weight (Anglo-American): 10.2 lb. Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Non-autoclavable Closure Size: 38 to 430 mm Diameter (metric): 59 mm Height (Anglo-American): 5.75 in. Opening: Narrow Width (anglo): 12.12 in. Unit Size: Case of 48

Additional Information

No.	CD-3637
Manufacturer (Brand)	Nalgene

