



Nalgene® square PETG media bottles

Product Images

Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 60 mL Capacity (Anglo): 2 oz. Diameter (anglo): 1.6 in. Inner diameter (anglo) neck: 0.71 in. Inner diameter (metric) neck: 18 mm For use with (application): Sampling, Preparation, Containment, Prolonged Storage Number per carton: 72 Number per package: 6 Temperature range (metric): gt;40deg;C Weight (anglo): 4.6 lb Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Not autoclavable Closure Size: 24 mm diameter (metric): 41 mm Height (anglo): 3.22 in. Height (metric): 82 mm Opening: Narrow Width (Anglo): 7.6 in. Unit Size: Case of 72

Additional Information

No.	CD-3635
Manufacturer (Brand)	Nalgene

