



Nalgene® square PETG media bottles

Product Images

Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 1 L Capacity (anglo): 32 oz. Diameter (anglo): 3.7 in. Inner diameter (anglo) neck: 1.10 in. Inner diameter (metric) neck: 28 mm For use with (application): Sampling, Preparation, Containment, Prolonged Storage Height (metric) with closure: 219 mm Number per carton: 35 Number per package: 5 Temperature range (metric): gt;40deg;C Weight (Anglo-American): 17.5 lb. Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Non-autoclavable Closure Size: 38 to 430 mm Height (Anglo-American): 8.67 in. Opening: Narrow Width (Anglo-American): 20.5 in. Unit Size: Case of 35

Additional Information

No.	CD-3639
Manufacturer (Brand)	Nalgene

