



Nalgene® square PETG media bottles

Product Images

Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 500 mL Capacity (Anglo): 16 oz. Diameter (anglo): 2.9 in. Inner diameter (Anglo-American) neck: 1.10 in. Inner diameter (metric) neck: 28 mm For use with (application): Sampling, Preparation, Containment, Prolonged Storage Height (metric) with closure: 177 mm Number per carton: 70 Number per package: 5 Temperature range (metric): gt;40deg;C Weight (Anglo-American): 7.6 lb. Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Non-autoclavable Closure Size: 38 to 430 mm Height (Anglo-American): 6.97 in. Aperture: Narrow Width (Anglo-American): 15.6 in. Unit Size: Case of 70

Additional Information

No.	CD-3638
Manufacturer (Brand)	Nalgene

