



# Nalgene® square PETG media bottles

## Product Images

## Description

Nalgene square PETG media bottles Thermo Scientific Nalgene square PETG media bottles are certified sterile and non-pyrogenic. Capacity (metric): 2000 mL Capacity (anglo): 64 oz. Diameter (anglo): 4.6 in. Inner diameter (Anglo-American) neck: 1.50 in. Inner diameter (metric) neck: 39 mm For use with (application): Sampling, Preparation, Containment, Prolonged Storage Height (metric) with closure: 271 mm Number per carton: 20 Number per package: 5 Temperature range (metric): gt;40deg;C Weight (Anglo-American): 14.5 lb. Certification/Compliance: USP Class VI Compliant, Certified Radiation Stabilized (SAL 10sup-6/sup) Autoclavable: Non-autoclavable Closure Size: 53B Diameter (metric): 116 mm Height (anglo): 10.67 in. Height (metric): 271 mm Opening: Narrow Width (Anglo-American): 19.75 in. Unit Size: Case of 20

## Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-3640 |
| Manufacturer (Brand) | Nalgene |

