



Nalgene® Autoclavable Syringe Bottles, PPCO

Product Images



Description

Nalgene Autoclavable Syringe Bottles, PPCO Thermo Scientific Nalgene autoclavable polypropylene copolymer syringe bottles are suitable for dispensing sterile water or solutions when the risk of cross-contamination must be minimized. Capacity: 1000 ml Capacity (Anglo-American): 32 ounces Includes: With closure/riser holder and riser Material: PPCO Number Per Case: 12 Number Per Package: 4 Type: Wash Bottle Closure Included: Yes Autoclavable: Autoclavable Closure Size: 38 to 430 mm Closure Material: Polypropylene Closure Type: Cap Graduated: No Height (Anglo-American): 10.2 in. Opening: Narrow Shape: Cylindrical Dispensing through Closure: Yes Unit Size: Case of 12

Additional Information

No.	CD-0600
Manufacturer (Brand)	Nalgene
VGKL number	228724052

