



Nalgene® Autoclavable Syringe Bottles, PPCO

Product Images



Description

Nalgene Autoclavable Syringe Bottles, PPCO Thermo Scientific Nalgene autoclavable polypropylene copolymer syringe bottles are suitable for dispensing sterile water or solutions when the risk of cross-contamination must be minimized. Capacity: 16 oz. (500 ml) capacity (Anglo-American): 16 oz. Color: Translucent Includes: With closure/riser holder and riser Material: PPCO Number Per Carton: 24 Number Per Package: 6 Type: Rinsing Bottle Closure Included: Yes Autoclavable: Autoclavable Closure Size: 28 mm Closure Material: Polypropylene Closure Type: Cap Graduated: No Height (Anglo-American): 8.4 in. Opening: Narrow Shape: Cylindrical Dispensing through closure: Yes Unit Size: Case of 24

Additional Information

No.	CD-0599
Manufacturer (Brand)	Nalgene
VGKL number	228724051

