



Nalgene® Teflon™ FEP syringe bottles

Product Images



Description

Nalgene™ Teflon™ FEP Squirt Bottles Ensure contamination-free dispensing with chemically resistant Thermo Scientific Nalgene Teflon FEP flush bottles molded entirely from fluoropolymers. Capacity: 4 oz., 125 mL capacity (metric): 125 mL capacity (anglo): 4 oz. Color: Translucent Type: Rinse Bottle Closure Included: Yes Autoclavable: Autoclavable Closure Size: 24 mm Closure Material: ETFE Closure Type: Screw Cap Graduated: No Height (Anglo): 6.2 in. Opening: Narrow Shape: Cylindrical Dispensing through closure: Yes Unit Size: Case of 4

Additional Information

No.	CD-0595
Manufacturer (Brand)	Nalgene
VGKL number	228724031

