



Nalgene® Unitary™ LDPE Squirt Bottles

Product Images



Description

Nalgene™ Unitary™ LDPE Syringe Bottles Dispense contents without tipping or shaking with Thermo Scientific™ Nalgene™ Unitary™ LDPE Syringe Bottles, which are faster and easier to fill with their wide neck. Capacity: 750 ml Capacity (Anglo): 24 ounces Description: LDPE bottle/tube neck Quantity per case: 24 Quantity per package: 4 Type: Wash Bottle Closure Included: Yes Closure Size: 38 mm Closure Material: Polypropylene Closure Type: Screw Cap Graduated: No Height (Anglo): 8.05 in. Opening: Wide Shape: Cylindrical Dispensing through Closure: No Unit Size: Case of 24

Additional Information

No.	CD-0593
Manufacturer (Brand)	Nalgene
VGKL number	228724024

