



Nalgene® Dropper Bottles

Product Images



Description

Nalgene Dropper Bottles For accurate, reliable and repeatable dispensing of various reagents. Made of opaque LDPE, these dropper bottles are ideal for light-sensitive solutions. An outstanding alternative to pipettors. Capacity (metric): 15 ml Capacity (anglo): 0.5 oz. Description: LDPE dropper bottle Contains: Dispensing tip and closures Quantity per carton: 25 Quantity per package: 25 Type: Dropper bottle Closure Color: White Certification/Compliance: ISO 13485:2003 Closure Included: Yes Closure Size: 15 to 415 mm Closure Type: drip control tip Dropper Capacity: 15 ml Drop Size: 40 µl Shape: Round Dispensing through Closure: Yes Unit Size: Case of 25

Additional Information

No.	CD-0702
Manufacturer (Brand)	Nalgene
VGKL number	228719050

