



Nalgene® Dropper Bottles

Product Images



Description

Nalgene Dropper Bottles For accurate, reliable and repeatable dispensing of various reagents. Made of opaque LDPE, these dropper bottles are ideal for light-sensitive solutions. An outstanding alternative to pipettors. Capacity (metric): 15 ml Capacity (anglo): 0.5 oz. Description: LDPE dropper bottle Contains: Dispensing tip and caps Quantity per case: 25 Quantity per package: 25 Type: Dropper bottle Cap color: Assorted colors Closure Included: Yes Closure Size: 15 to 415 mm Closure Type: Drip Control Tip Drop Size: 40 µl For Use With (Devices): Nalgene closures Shape: Round Dispensing through Closure: Yes Unit Size: Case of 25

Additional Information

No.	CD-0707
Manufacturer (Brand)	Nalgene
VGKL number	228739050

