




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Matrix ScrewTop Tri-coded Tubes

€486.00
plus VAT &
Shipping

Product Images

Description

Matrix ScrewTop Tri-coded Tubes Securely store precious samples at low temperatures, including vapor phase Liquid Nitrogen with Thermo Scientific Matrix 500l and 1.0mL ScrewTop Tubes. Capacity (Anglo-American): 0.0338 fl oz. DNase-free: Yes Type: Screw-Top Tubes Description: V-Bottom ScrewTop Tubes with white patch and linear barcode Cap Color: Colorless Capacity (metric): 1.0 mL Bottom shape: V-bottom Pack configuration: Bulk Well shape: V-bottom Barcode: Non-proprietary, 12x12 Data-matrix with ECC200 Built-in Error Correction Unit Size: Each

Additional Information

No.	TA-0044
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

