



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Nunc™™ 2.0 ml universal tubes with internal thread

€633.60
plus VAT &
Shipping

Product Images



Description

Nunc™ 2.0 ml Universal Tubes with Internal Thread Thermo Scientific Nunc 2.0 ml Universal Tubes with Internal Thread are available with uncoded 2D barcodes, without barcodes, or both linear and 2D barcodes to meet all requirements. Now available in a hinged cap rack with barcode. Cap color: Natural Description: Nunc Universal Tubes with 2D code Capacity (metric): 2.0 mL Material: plastic Sterility: sterile SAL 10sup-6/sup Sterility: sterile SAL 10sup-6/sup Type: universal tubes Coding: 2D on bottom of tube Vial type: internal thread, flat bottom Color: Natural Colors Packaging: bulk Quantity per box: 480 For use with (devices): Biobank, Cryostorage Contains: White field and graduations Quantity: 480 Unit Size: Case of 480

Additional Information

No.	FE-0755
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

