



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Nunc™™ 2.0 ml universal tubes with internal thread

€173.00
plus VAT &
Shipping

Product Images



Description

Nunc™ 2.0 ml Universal Tubes with Internal Thread Thermo Scientific Nunc 2.0 ml Universal Tubes with Internal Thread are available with uncoded 2D barcodes, without barcodes, or both linear and 2D barcodes to meet all requirements. Now available in a hinged lid rack with barcode. Description: Nunc latch rack Capacity (metric): 2 ml (tubes) Material: Plastic Type: Universal Tubes, Universal Hinged Lid Rack Coding: None Vial Type: Threaded Universal Tubes Color: Beige with Blue Latches Packaging: 10 empty racks Number per Package: 1 Number per Case: 10 Quantity: 10 Unit Size: Case of 10

Additional Information

No.	FE-0600
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

