



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Open Top Matrix Storage Tubes with 2D Barcode in 384 Format

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Open top matrix storage tubes with 2D barcode in 384 format Maximize storage capacity with the higher density 384 tube storage format. Capacity (anglo): 0.0034 fl oz. Package Quantity: 20 carton (racks of 384 tubes) Type: Open Top Type: 2D Barcode Package Configuration: 20 racks of 384 tubes / cs . Well Shape: Flat Barcode: 2D BarCoded Product Line: Matrixtrade; Unit Size: Case of 7680

Additional Information

No.	TA-0124
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

