



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific tubes with linear barcode

€1,800.00
plus VAT &
Shipping

Product Images



Description

Linear Barcode Tubes Simplify cryogenic tube tracking with linear barcodes printed directly on the side of the tubes. Single-handed operation: Yes Capacity (Anglo-American): 0.06086 oz. Cap type: Screw cap DNase-free: Yes Package quantity: 50 packs, 1,800 boxes Type: Linear barcode Certification/Compliance: Nunc tubes are compliant with IATA for shipping. CE Marked for storage of diagnostic samples in the EU. Description: Nunc Barcoded (inkjet interleave 2 of 5) CryoTube Vial Autoclavable: Not autoclavable Thread type: internal Length (metric): 48. mm Quantity per carton: 1800 Quantity per pack: 50 RNase-free: yes Barcode: Yes Product Line: Nunctrade; Unit Size: Case of 1800

Additional Information

No.	FE-0700
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

