



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand® cuvette, 4 optical windows, macro, nest  
number sorted, PS, minimum filling volume 2.5  
ml**

**€71.73** As low as  
**€59.83**  
**plus VAT &  
Shipping**

## Product Images



## Description

---

Disposable cuvettes 4 optical windows, macro, PS or UV transparent, nest number sorted, minimum filling volume 2.5 ml. Suitable for fluorescence spectroscopy. The UV version can be used from a wavelength of 230 nm and shows no autofluorescence. Quality features: Optically perfect transmission range. Production under controlled room conditions and fully automated packaging. Sorted by nest numbers to ensure the lowest possible deviations in absorbance values. The UV-transparent cells replace sensitive and expensive glass or quartz glass cells in many areas.

## Additional Information

---

No.	BN-2093
Manufacturer (Brand)	Brand
VGKL number	104975035
Manufacturer EAN	4033378488398
Transport temperature	Room temperature

