



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand® cuvette, 4 optical windows, macro,
nest-number-sorted, PS, minimum filling
volume 2.5 ml, 759030**

€15.43 As low as
€12.89
**plus VAT &
Shipping**

Product Images



Description

Disposable cuvettes 4 optical windows, macro, PS or UV transparent, nest number sorted, minimum filling volume 2.5 ml. Suitable for fluorescence spectroscopy. The UV version can be used from a wavelength of 230 nm and shows no autofluorescence Quality features: Optically perfect transmission range Production under controlled room conditions and fully automated packaging Sorted by nest numbers to ensure the lowest possible deviations in absorbance values The UV-transparent cells replace sensitive and expensive glass or quartz glass cells in many areas

Additional Information

No.	BN-2092
Manufacturer (Brand)	Brand
VGKL number	104975030
Manufacturer EAN	4033378488428
Transport temperature	Room temperature

