



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Liquid Filter Set UV101H

€1,091.00
plus VAT &
Shipping

Product Images

Description

Liquid filter set type 667-UV101H for checking scattered light according to USP 857, filter set consisting of 2 liquid filters: Sodium iodide in ultrapure water with layer thickness 10 mm (UV10) and reference filter sodium iodide in ultrapure water with layer thickness 5 mm (UV10H), certified at wavelength: 259 nm (cut-off), in high-quality storage box with DAkKS calibration certificate Attention: article available until 30.06.2022. Further information will follow.

Additional Information

No.	KI-0354
Manufacturer (Brand)	Hellma
Transport temperature	Room temperature

