



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma liquid filter set 667-UV100H

€1,091.00
plus VAT &
Shipping

Product Images

Description

Liquid filter set type 667-UV100H for checking scattered light according to USP 857, filter set consisting of 2 liquid filters: potassium chloride in ultrapure water with layer thickness 10 mm (UVI) and reference filter potassium chloride in ultrapure water with layer thickness 5 mm (UVIH), certified at wavelength: 200 nm (cut-off), in high-quality storage box with DAkkS calibration certificate Attention: Article available until 30.06.2022. Further information will follow.

Additional Information

| | |
|-----------------------|------------------|
| No. | KI-0352 |
| Manufacturer (Brand) | Hellma |
| Transport temperature | Room temperature |

