



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma glass filter set 666-S300

€1,150.00
plus VAT &
Shipping

Product Images

Description

Glass filter set type 666-S300 for checking the photometric accuracy consisting of 3 neutral glass filters (F390, F301, F303) with the absorptions: 0.04 / 2.5 / 3.0, certified at the wavelengths: 440; 465; 546.1; 590; 635 nm, in high-quality storage box with DAkkS calibration certificate Attention: Article available until 30.06.2022. Further information will follow.

Additional Information

No.	KI-0316
Manufacturer (Brand)	Hellma
Transport temperature	Room temperature

