



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Neutral Glass Filter 666-F303 (3.0 Abs)

€410.00
plus VAT &
Shipping

Product Images

Description

Neutral density glass filter type 666-F303 for checking the photometric accuracy (3.0 Abs) certified at the wavelengths: 440; 465; 546.1; 590; 635 nm, in high-quality storage box with DAkKS calibration certificate Attention: article available until 30.06.2022. Further information will follow.

Additional Information

No.	KI-0683
Manufacturer (Brand)	Hellma
Transport temperature	Room temperature

