



Thermo Scientific 96-Well Microtiter Microplates

Product Images



Description

96-Well Microtiter Microplates Use these 96-well microtiter plates for microvolume applications such as storage, sample or reagent transfer, serial dilution, hemagglutination, blood typing, and complement fixation. Bottom: C Format: Standard Lid: Without lid Array: 8 x 12 Packaging: 50 Cartons, 10 Pkg. Product Line: Microtitertrade; Product Type: Microtiter Plate Volume (metric) Well: 330 µl Number of Wells: 96 Surface Treatment: Untreated Number per Carton: 50 Number per Package: 10 Quantity: 50 pcs./Pack Unit Size: Case of 50

Additional Information

No.	FE-0950
Manufacturer (Brand)	Thermo Scientific

