




Thermo Scientific Nunc™ microplate coverslips

Product Images



Description

Nunc Microplate Lids Protect samples from contamination and evaporation during assay processing, incubation or storage with Thermo Scientific™ Nunc™ microplate lids. Dimensions (L x W x H): 5 x 3.35 x 0.18 in. (127 x 85 x 4.6 mm) Type: lowbot lid Number of wells: 384/1536 Length (anglo): 5 in. Length (metric): 127 mm Height (anglo): 0.18 in. Quantity per carton: 100 Quantity per package: 25 Width (anglo): 3.35 in. Plate Compatibility: 384-well plates, 1536-well plates Condensation rings: Yes Evaporation barrier: Yes Chamfered corners: No For use with: 384- and 1536-well plates Unit Size: Case of 100

Additional Information

No.	FE-0462
Manufacturer (Brand)	Thermo Scientific

