



Thermo Scientific Nunc™ NucleoLink/TopYield Frame

Product Images



Description

Nunc NucleoLink/TopYield Frames The unique design of these frames is ideal for liquid phase PCR, immuno-PCR and hybridization assays, providing error-free flexibility, accurate reading and efficient washing. Base: V Format: Standard Lid: Without Lid Description: Nunc NucleoLink/TopYield Frames Array: 8 x 12 Packaging: 72 boxes Product Type: Frames Material: ABS plastic Number of Wells: 96 Well Shape: V-shaped For use with (application): Liquid phase PCR or immuno-PCR Max. Heat resistance: 90 °C Number per carton: 72 Number per package: 6 Plate compatibility: 96-well plates Quantity: 72 Length (metric): 128 mm Width (metric): 86 mm Unit Size: Case of 72

Additional Information

No.	FE-0425
Manufacturer (Brand)	Thermo Scientific

