



Thermo Scientific Nunc™™ MicroWell™ MiniTrays

Product Images

Description

Nunc™ MicroWell™ MiniTrays More space in freezers and incubators with Nunc MicroWell MiniTrays, which can be used with inverted microscopes thanks to their optical transmission. Closure included: yes Closure material: Polystyrene Description: MicroWell MiniTray for serological applications Shape: Rectangular Surface treatment: Cell culture treated Type: Cell Culture Dish Number of wells: 72 For use with (application): For serological applications Contains: Lid Number per box: 100 Number per package: 10 Height (metric): 59 mm Unit Size: Case of 100

Additional Information

No.	FE-0805
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900438733

