



Nalgene® Polycarbonate Centrifuge Bottles

Product Images



Description

Nalgene™ Polycarbonate Centrifuge Bottles Use these polycarbonate centrifuge bottles for harvesting, pelleting, and protein purification of large cell volumes. Capacity (metric): 500 mL Capacity (anglo): 16.9 oz. Description: polycarbonate centrifuge bottles Number per carton: 24 Number per package: 4 Type: centrifuge bottle Nominal capacity: 500 mL Closure Included: Yes Autoclavable: autoclavable Closure material: polypropylene Closure type: screw cap For use with (devices): Fiberlitrade; rotors Graduated: No Max. Force: 13,700 x G Shape: Round Dispensing Through Closure: No Unit Size: Case of 24

Additional Information

No.	CD-0934
Manufacturer (Brand)	Nalgene
VGKL number	238731223

