



Nalgene® Polycarbonate Centrifuge Bottles

Product Images



Description

Nalgene™ Polycarbonate Centrifuge Bottles Use these polycarbonate centrifuge bottles for harvesting, pelleting, and protein purification of large cell volumes. Capacity (metric): 1000 mL Capacity (anglo): 33.81 oz. Description: polycarbonate centrifuge bottles Number per carton: 16 Number per pack: 4 Type: centrifuge bottle Nominal capacity: 1000 mL Autoclavable: autoclavable Closure material: polypropylene Closure type: screw cap For use with (devices): Fiberlitrade; rotors Graduated: No Max. Force: 7100 x g Shape: Round Dispensing through the cap: No Unit Size: Case of 16

Additional Information

No.	CD-0935
Manufacturer (Brand)	Nalgene
VGKL number	238731224

