



Nalgene® Teflon FEP Centrifuge Bottles

Product Images

Description

Nalgene Teflon FEP Centrifuge Bottles Thermo Scientific Nalgene centrifuge bottles are ideal for lipid and phenol extractions, digestions and precipitations at low speeds. Capacity (metric): 250 ml Capacity (anglo): 8.45 oz. Description: Teflon FEP bottle with Tefzel ETFE closure Number per carton: 6 Number per package: 1 Type: Centrifuge bottle Nominal capacity: 250 ml Closure material: Tefzel ETFE Closure type: Screw cap Graduated: No Max. Force: 4000 x G Shape: Round Dispensing Through Cap: No Unit Size: Case of 6

Additional Information

No.	CD-0940
Manufacturer (Brand)	Nalgene

