




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Wide Neck Centrifuge Bottles, HDPE

€313.92
plus VAT &
Shipping

Product Images

Description

Nalgene wide neck centrifuge bottles, HDPE Thermo Scientific Nalgene centrifuge bottles are economical products. Capacity (metric): 250 mL Capacity (anglo): 8.45 oz. Description: high density polyethylene bottle with PP cap Material: HDPE (high density polyethylene) Length (metric): 132 mm Quantity per carton: 36 Quantity per package: 6 Type: bottle Outer diameter (metric): 61 mm Nominal capacity: 250 ml Closure material: polypropylene Closure type: screw cap Graduated: no Max. Force: 8000 × g Opening: Wide Shape: Round Dispensing through cap: No Unit Size: Case of 36

Additional Information

No.	CD-0873
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

