



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® centrifuge bottles made of PPCO

€52.40 As low
as
€42.24
plus VAT &
Shipping

Product Images



Description

Nalgene PPCO centrifuge bottles Perform biopharmaceutical aseptic liquid transfers or sampling with the three-sided flange of Thermo Scientific Nalgene PPCO centrifuge bottles. Capacity (metric): 250 mL Capacity (anglo): 8.45 oz. Color: Translucent Material: Polypropylene copolymer (PPCO) Quantity per case: 36 Quantity per package: 4 Type: Centrifuge bottle Nominal capacity: 250 mL Autoclavable: Autoclavable Closure material: Polypropylene (PP) Closure type: Screw cap Graduated: No Max. Force: 13,200 x g (up to 27,500 x g with seal closure unit) Shape: Round Dispensing through cap: No Unit Size: Case of 36

Additional Information

No.	CD-0778
Manufacturer (Brand)	Nalgene
VGKL number	238731201
Transport temperature	Room temperature

