




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene® centrifuge bottles made of PPCO

**€60.00** As low  
as  
**€48.32**  
plus VAT &  
Shipping

### Product Images



## Description

---

Nalgene PPCO centrifuge bottles Perform biopharmaceutical aseptic liquid transfers or sampling with the three-sided flange of Thermo Scientific Nalgene PPCO centrifuge bottles. Capacity (metric): 500 mL Capacity (Anglo): 16.9 oz. Color: Translucent Material: Polypropylene copolymer (PPCO) Quantity per case: 24 Quantity per package: 4 Type: Centrifuge bottle Nominal capacity: 500 mL Autoclavable: Autoclavable Closure material: Polypropylene (PP) Closure type: Seal closure Graduated: No Max. Force: 13,700 × g Shape: Round Dispensing Through Closure: No Unit Size: Case of 24

## Additional Information

---

No.	CD-0861
Manufacturer (Brand)	Nalgene
VGKL number	238731203
Transport temperature	Room temperature

