




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Tubes and Stoppers for Immunology

€216.00
plus VAT &
Shipping

Product Images



Description

Tubes and plugs for immunology With Thermo Scientific™ Immuno Tubes and Plugs, you can reproduce your results over and over again. They are also ideal for ELISA, solid phase RIA and liquid phase RIA assays. Material: Polyethylene Shape: Square Type: Immuno Tubes and Stoppers Description: liquid phase polyethylene tubes Finish: MiniSorp Total Volume: 12 mL Number Per Package: 150 Number Per Carton: 1200 Volume (metric): 12 mL For use with (devices): Immunology Tubes Unit Size: Case of 1200

Additional Information

No.	FE-0850
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900468608
Transport temperature	Room temperature

