




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Tubes and Stoppers for Immunology

€330.00
plus VAT &
Shipping

Product Images



Description

Tubes and plugs for immunology With Thermo Scientific™ Immuno Tubes and Plugs, you can reproduce your results over and over again. They are also ideal for ELISA, solid phase RIA and liquid phase RIA assays. Material: Polyethylene Shape: Square Type: Immuno Tubes and Stoppers Description: liquid phase polyethylene Tubes Finish: MiniSorp Total Volume: 5 mL Number per Package: 100 Number per Carton: 3000 Volume (metric): 5 mL For Use with (Devices): Immunology Tubes Unit Size: Case of 3000

Additional Information

No.	FE-0811
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

