



## Hirschmann® Silicoses, pink T28

### Product Images

---



## Description

---

Silicosen is sterile stoppers produced by a special process, which have regular pores and are therefore ideal for the preparation and sterilization of culture media. The uniform pores and good air permeability make these stoppers ideal for culturing aerobic microorganisms. These designs are chemical and temperature resistant as well as liquid repellent and reusable.

## Additional Information

---

|                      |            |
|----------------------|------------|
| No.                  | NH-1071    |
| Manufacturer (Brand) | Hirschmann |
| VGKL number          | 103705528  |

