



Nalgene® microtiter tubes and racks

Product Images

Description

Microtiter Tubes and Racks These polyethylene plug bars for use with 1.2mL MCT microtiter racks are designed to tightly seal microtiter tubes for proper storage. Material: polyethylene (stopper) Type: microtiter Certification/Compliance: Certified free of RNase, DNase, DNA and pyrogens Color: Natural Description: 1.2mL MicroTiter Tube Robotic Rack, Nonsterile Volume (metric): 1.2 mL Quantity: 96 tubes/tray 10 trays/pack, nonsterile Product Line: Thermo Scientifictrade; Unit Size: Case of 4800

Additional Information

No.	CD-3216
Manufacturer (Brand)	Nalgene

