



Nalgene® microtiter tubes and racks

Product Images

Description

Microtiter tubes and racks These polyethylene plug bars for use with 1.2 ml MCT microtiter racks are used to tightly seal microtiter tubes for proper storage. Material: polyethylene (stopper) Type: microtiter Certification/Compliance: certified free of RNase, DNase, DNA and pyrogens Color: natural Description: 1.2 ml microtiter tubes, bulk pack, non-sterile Volume (metric): 1.2 ml Quantity: 1000 tubes/unit, 10 units/case Product Line: Thermo Scientifictrade; Unit Size: Case of 10,000

Additional Information

No.	CD-3215
Manufacturer (Brand)	Nalgene

