



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® PPCO Centrifuge Tubes with Round Bottom for

€392.80
plus VAT &
Shipping

Product Images



Description

Nalgene PPCO round-bottom centrifuge tubes for high-speed applications Thermo Scientific™ Nalgene™ PPCO round-bottom centrifuge tubes for high-speed applications have uniform wall thickness and diameter, fit better in rotors, and perform better. Capacity (metric): 100 mL Quantity per box: 100 Quantity per pack: 10 Autoclavable: yes Graduations: No Shape: Round Type: Centrifuge Tube Max. Force: 50,000 x g Writing Area: No Unit Size: Case of 40

Additional Information

No.	CD-0744
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

