



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® PPCO Centrifuge Tubes with Round Bottom for

€375.00
plus VAT &
Shipping

Product Images



Description

Nalgene PPCO round-bottom centrifuge tubes for high-speed applications Thermo Scientific™ Nalgene™ PPCO round-bottom centrifuge tubes for high-speed applications have uniform wall thickness and diameter, fit better in rotors, and perform better. Material: polypropylene copolymer (PPCO) Quantity per box: 100 Quantity per pack: 10 Autoclavable: Yes Graduations: No Shape: Round Type: Centrifuge Tube Max. Force: 50,000 x g Writing Area: No Unit Size: Case of 100

Additional Information

| | |
|-----------------------|------------------|
| No. | CD-0739 |
| Manufacturer (Brand) | Nalgene |
| Transport temperature | Room temperature |

