



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® centrifuge tubes made of polycarbonate and with round bottom for

**€447.00
plus VAT &
Shipping**

Product Images



Description

Nalgene polycarbonate and round-bottom centrifuge tubes for high-speed applications The exceptional mechanical strength and high optical clarity of these polycarbonate round-bottom centrifuge tubes allow you to successfully centrifuge samples multiple times. Capacity (metric): 16 ml Quantity per box: 100 Quantity per pack: 10 Autoclavable: yes Graduations: No Shape: Round Type: Centrifuge Tube Max. Force: 50,000 x g Writing Area: No Unit Size: Case of 100

Additional Information

No.	CD-0757
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

