



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Nalgene® centrifuge tubes made of polycarbonate and with round bottom for**

**€557.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

Nalgene polycarbonate and round-bottom centrifuge tubes for high-speed applications Thanks to the exceptional mechanical strength and high optical clarity of these polycarbonate round-bottom centrifuge tubes, you can successfully centrifuge samples multiple times. Capacity (metric): 50 ml Number per box: 100 Number per pack: 10 Autoclavable: Yes Graduations: No Shape: Round Type: Centrifuge Tube Max. Force: 50,000 x G Length (metric): 106.0 mm Writing Area: No Unit Size: Case of 100

## Additional Information

---

No.	CD-0761
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

