



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton® 4 ml sample vial, brown borosilicate glass type I, black phenolic resin cap, 13-425, without hole

**€92.40
plus VAT &
Shipping**

Product Images

Description

4 ml sample vial, brown borosilicate glass type I, with mounted black phenolic cap 13-425 without hole, with 14B rubber insert (pack of 40 pcs.) DURAN code: 4027745

Additional Information

No.	WD-2147
Manufacturer (Brand)	Wheaton
Transport temperature	Room temperature

