



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton® 4ml Sample Vial Lab File Borosilicate glass type I, black, phenolic resin cap, PTFE-coated

€93.60
plus VAT &
Shipping

Product Images

Description

4 ml sample vials in Lab File®, clear borosilicate glass type I, with mounted black phenolic cap 13-425 without hole, PTFE coated 14B rubber insert (pack of 40 pcs.) DURAN code: 4027623

Additional Information

| | |
|-----------------------|------------------|
| No. | WD-2137 |
| Manufacturer (Brand) | Wheaton |
| Transport temperature | Room temperature |

