



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® threaded tube for attachment, GL 32

€107.50
plus VAT &
Shipping

Product Images



Description

DURAN® threaded tube for attachment, GL 32 This threaded tube has a screw thread and can be incorporated into glass components by a further processing step. Special properties and advantages: With DIN thread Application: Attachment in glass apparatus

Additional Information

No.	SD-1694
Manufacturer (Brand)	DURAN
VGKL number	132483901
Transport temperature	Room temperature

