



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Vacuum distillation pusher, straight, with unscrewable plastic olive, NS 29/32

€60.52
plus VAT &
Shipping

Product Images



Description

DURAN® Vacuum distillation spur, straight, NS 29/32 One part of a distillation setup is the distillation spur, which is connected downstream at the end of the cooler and feeds the distillate into the receiving flask. If vacuum distillation is performed, the vacuum pump is ideally connected to the distillation head. Special features and advantages: Tight closure due to standard ground joint Straight With two NS and unscrewable plastic olive Burnt-in and therefore permanently durable white print Application: Distillation

Additional Information

No.	SD-1140
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

