



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Intensive Cooler, with 2 unscrewable plastic olives, NS 24/29, jacket length 400 mm

€141.70
plus VAT &
Shipping

Product Images



Description

DURAN® Intensive Cooler, NS 24/29, 400 mm Combining the cooling properties of the Liebig cooler with those of the Dimroth cooler produces the intensive cooler. This is characterized by its large cooling surface (cooling coil and outer jacket) and is thus suitable for condensing vapors from reaction mixtures with low boiling points (distillation). Special features and advantages: This type of cooler has a very large cooling surface due to a cooling coil and a double jacket and is therefore particularly suitable for work with low-boiling media Two unscrewable plastic olives for inflow and outflow of the cooling medium Tight closure due to standard ground joint Application: Reflux condenser for condensation and return of the (solvent) vapors to the reaction mixture, distillation, All coolers are special sizes and not according to DIN.

Additional Information

No.	SD-1291
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

