



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® ball cooler, with 2 unscrewable plastic olives, NS 24/29, jacket length 250 mm

**€90.45**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® ball cooler (Allihn cooler), NS 24/29, 250 mm The ball cooler is similar in design to the Liebig cooler. Its internal riser tube is constricted at equal intervals. This creates the typical spherical shape that gives the cooler its name. Due to this shape, the cooling surface is enlarged compared to the Liebig cooler and produces a stronger steam turbulence. Special features and advantages: Compared to the Liebig cooler, the spherical cooler has a larger cooling surface and thus a higher cooling capacity Two unscrewable plastic olives for inflow and outflow of the cooling medium Tight closure due to standard ground joint Application: reflux condenser for condensation and return of the (solvent) vapors to the reaction mixture, distillation, product and reflux condenser, coolers of length 250 mm (with NS 24/29 and NS 29/32) and 400 mm (with NS 24/29) are special sizes and not according to DIN.

## Additional Information

---

No.	SD-1279
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

