



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Liebig cooler, with 2 unscrewable plastic olives, NS 29/32, sheath length 400 mm

**€71.16**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

DURAN® Liebig condenser (West condenser), NS 29/32, 400 mm The Liebig condenser is one of the best known types of condenser and is used for the distillation of reaction mixtures. The cooler is a tube with standard ground joints at both ends and a glass jacket. The cooling water passed through the jacket in counterflow principle to increase the cooling capacity. Special features and advantages: Small cooling surface and therefore relatively low cooling capacity Two unscrewable plastic olives for inflow and outflow of the cooling medium Tight seal due to standard ground joint Application: product cooler for separation of the distillate, distillation, cooling in countercurrent principle, coolers of length 250 mm (with NS 14/23, NS 24/29 and NS 29/32) and 400 mm (with NS 24/29) are special sizes and not according to DIN.

## Additional Information

---

No.	SD-1278
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

