



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**DURAN® snake cooler, with 2 unscrewable plastic olives, NS 29/32, jacket length 300 mm**

**€114.39  
plus VAT &  
Shipping**

## Product Images



## Description

---

DURAN® serpentine cooler, NS 29/32, 300 mm Typical of the serpentine cooler is its coiled riser tube, which is located in a jacket and is thus surrounded by the cooling medium (contrary to the Dimroth cooler). Due to this condensate flow, the condensation capacity is very low although the cooling distance is very long. This type of cooler is therefore usually used as a downstream additional cooler and not as an independent distillation cooler. Due to the better drainage of the condensate in the coiled riser tube, the cooler design is required in vertical position. Special features and advantages: Due to the long cooling distance mostly used as a downstream additional cooler Two unscrewable plastic olives for inflow and outflow of the cooling medium Tight closure due to standard ground joint Application: Cooling of a condensate, product cooler for separation of the distillate, distillation

## Additional Information

---

No.	SD-1284
Manufacturer (Brand)	DURAN
Transport temperature	Room temperature

