



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® bells for negative pressure chambers

€264.00
plus VAT &
Shipping

Product Images

Description

Nalgene Bells for Negative Pressure Chambers Thermo Scientific Nalgene polyetherimide replacement bells for negative pressure chambers are chemical resistant and have an amber color. Closure Type: Brim Cap For use with (devices): Nalgenetrade; Vacuum Chambers Capacity (Anglo-American): 2.625 gal. Capacity (metric): 8.3L Diameter (Anglo): 8.75 in. Height (anglo): 10 in. Material: Polyetherimide Type: Replacement Vacuum Jars ; 2.13 gal . (8.3L) Unit Size: Each

Additional Information

No.	CD-2178
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

